



### **Product description**

|                     |  |
|---------------------|--|
| Order number        | 72.985.002   |
| Product description | PCR strip of 8, up to 200 µl, High Profile, PCR Performance Tested, transparent, material: PP, without cap, suitable for qPCR, 120 piece(s)/minigrip bag |

### **Product characteristics**

|                            |                |
|----------------------------|----------------|
| Graduation                 | no             |
| Suitable for qPCR          | yes            |
| Base shape                 | conical        |
| High or low profile        | High Profile   |
| Number of wells            | 8              |
| Volume of work             | 10 µl - 0,2 ml |
| Centrifugation max (RCF)   | 8000 x g       |
| Centrifugation temperature | 22             |
| Centrifugation time        | 5 Min          |
| Autoclavable (max. °C)     | 121 °C         |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

|                   |         |
|-------------------|---------|
| Sample volume     | 0.2 ml  |
| Max. volume       | 200 µl  |
| Diameter          | 5.4 mm  |
| Width of product  | 8.1 mm  |
| Height of product | 20.7 mm |
| Length of product | 70.9 mm |

### **Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | transparent        |

### **Purity & certification**

|                  |                                |
|------------------|--------------------------------|
| Product category | In vitro diagnostic device, CE |
| CE certified     | CE                             |
| Purity standard  | PCR Performance Tested         |
| Batched          | yes                            |

### **Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 480           |
| Type of smallest subpackaging | minigrip bag  |
| Piece(s) / inner box          | 120           |
| Piece (s) / outer case        | 480           |
| Piece(s) / pallet             | 155520        |
| Depth of case                 | 268 mm        |
| Width of case                 | 191 mm        |
| Height of case                | 100 mm        |
| Case volume                   | 0.0051 cbm    |
| Weight of product             | 0.001 kg      |
| Weight of case                | 0.63 kg       |
| EAN of inner box              | 4038917124027 |
| EAN case                      | 4038917019002 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

