



Product Specification
PCR strip of 4, 200 µl, Biosphere® plus,
transparent, PP, flat cap

Page 1



Product description

Order number	72.990
Product description	PCR strip of 4, up to 200 µl, High Profile, Biosphere® plus, transparent, material: PP, with anti-contamination shield, flat cap, cap attached, suitable for qPCR, 50 piece(s)/bag

Product characteristics

Graduation	no
Cap	flat cap
Closure type	snap-on lid
Suitable for qPCR	yes
Base shape	conical
High or low profile	High Profile
Number of wells	4
Volume of work	10 µl - 0,2 ml
Centrifugation max (RCF)	8000 x g
Centrifugation temperature	22
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
PCR strip of 4, 200 µl, Biosphere® plus,
transparent, PP, flat cap

Size

Sample volume	0.2 ml
Max. volume	200 µl
Diameter	5.35 mm
Width of product	10.6 mm
Height of product	20.65 mm
Length of product	37.2 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	400
Type of smallest subpackaging	bag
Piece(s) / inner box	200
Piece (s) / outer case	400
Piece(s) / pallet	57600
Depth of box	189 mm
Width of box	151 mm
Height of box	97 mm
Depth of case	323 mm
Width of case	213 mm
Height of case	145 mm
Case volume	0.01 cbm
Weight of product	0.0008 kg
Weight of case	0.63 kg
EAN of inner box	4038917084871

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
PCR strip of 4, 200 µl, Biosphere® plus,
transparent, PP, flat cap

Page 3

EAN case

4038917030045

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT