



Product description

| | |
|---------------------|---|
| Order number | 72.990.992 |
| Product description | PCR strip of 4, up to 200 µl, High Profile, PCR Performance Tested, red/green/blue/violet, material: PP, with anti-contamination shield, flat cap, cap attached, suitable for qPCR, 120 piece(s)/minigrip bag |

Product characteristics

| | |
|----------------------------|----------------|
| Graduation | no |
| Cap | flat cap |
| Closure type | push cap |
| Suitable for qPCR | yes |
| Base shape | conical |
| High or low profile | High Profile |
| Number of wells | 4 |
| Volume of work | 10 µl - 0,2 ml |
| Centrifugation max (RCF) | 8000 x g |
| Centrifugation temperature | 22 |
| Centrifugation time | 5 Min |
| Autoclavable (max. °C) | 121 °C |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

| | |
|-------------------|----------|
| Sample volume | 0.2 ml |
| Max. volume | 200 µl |
| Diameter | 5.35 mm |
| Width of product | 10.2 mm |
| Height of product | 20.65 mm |
| Length of product | 37.2 mm |

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | colour mix |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Purity standard | PCR Performance Tested |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 480 |
| Type of smallest subpackaging | minigrip bag |
| Piece(s) / inner box | 480 |
| Piece (s) / outer case | 480 |
| Piece(s) / pallet | 69120 |
| Depth of case | 323 mm |
| Width of case | 213 mm |
| Height of case | 145 mm |
| Case volume | 0.01 cbm |
| Weight of product | 0.0007 kg |
| Weight of case | 0.57 kg |
| EAN of inner box | 4038917124287 |
| EAN case | 4038917030113 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

