



### **Product description**

Order number	72.994.992
Product description	PCR strip of 8, High Profile, PCR Performance Tested, red/green/blue/violet/yellow, material: PP, round cap, cap attached, 120 piece(s)/bag

### **Product characteristics**

Graduation	no
Cap	round cap
Closure type	push cap
Base shape	conical
High or low profile	High Profile
Number of wells	8

### **Size**

Sample volume	0.2 ml
---------------	--------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material Polypropylene (PP)  
Colour of product colour mix

**Purity & certification**

Product category no medical device | IVD  
Purity standard PCR Performance Tested  
Batched yes

**Packaging**

Minimum order qty. 1200  
Type of smallest subpackaging bag  
Piece(s) / inner box 120  
Piece (s) / outer case 1200  
Piece(s) / pallet 1200  
Case volume 0.0299 cbm  
Weight of product 0.0031 kg  
Weight of case 4.3 kg  
EAN of inner box 4038917124331  
EAN case 4038917080743

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

