



Product description

Order number	95.1994
Product description	PCR film, free from DNase/RNase, material: PO, transparent, suitable for qPCR, self-adhesive, application area: -70 - 110 °C, strong adhesion, high chemical resistance, 100 piece(s)/bag

Product characteristics

Closure type	film
Light protection	no
Suitable for qPCR	yes
Minimum temperature	-70 °C
Maximum temperature	110 °C

Size

Width of product	79 mm
Length of product	141 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polyolefin (PO)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Purity standard free from DNase/RNase
Batched yes

Packaging

Minimum order qty. 100
Type of smallest subpackaging bag
Piece(s) / inner box 100
Piece (s) / outer case 1200
Piece(s) / pallet 288000
Case volume 0.0052 cbm
Weight of product 0.002 kg
Weight of case 2.72 kg
EAN of inner box 4038917088114
EAN case 4038917044301

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

