



Product description

Order number	95.1999
Product description	qPCR film, free from DNase/RNase, material: PO, highly transparent, suitable for qPCR, encapsulated adhesive, application area: -70 - 100 °C, heat sensitive adhesive, lowest evaporation rates, 100 piece(s)/bag

Product characteristics

Closure type	film
Light protection	no
Suitable for qPCR	yes
Minimum temperature	-70 °C
Maximum temperature	100 °C

Size

Width of product	77 mm
Length of product	140 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polyolefin (PO)
Colour of product	highly transparent

Purity & certification

Product category	no medical device IVD
Purity standard	free from DNase/RNase
Batched	yes

Packaging

Minimum order qty.	100
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1600
Piece(s) / pallet	192000
Depth of case	173 mm
Width of case	135 mm
Height of case	283 mm
Case volume	0.0066 cbm
Weight of product	0.0021 kg
Weight of case	3.92 kg
EAN of inner box	4038917303590
EAN case	4038917388238

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

