



### **Product description**

Order number	95.988
Product description	Base station, material: PC, transparent, (LxWxH): 128 x 85.5 x 30.4 mm, suitable for work tray 95.987.002 & 0.2 ml PCR tubes, 5 piece(s)/case

### **Product characteristics**

Suitable for	work tray 95.987.002 & 0.2 ml PCR tubes
--------------	---

### **Size**

Width of product	85.5 mm
Height of product	30.4 mm
Length of product	128 mm

### **Material & colours**

Product material	Polycarbonate (PC)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category                   no medical device | IVD  
Batched                               no

**Packaging**

Minimum order qty.               5  
Type of smallest                   case  
subpackaging  
Piece(s) / inner box               5  
Piece (s) / outer case           75  
Piece(s) / pallet                  3600  
Depth of case                    598 mm  
Width of case                     398 mm  
Height of case                    145 mm  
Case volume                      0.0345 cbm  
Weight of product                0.0302 kg  
Weight of case                    2.83 kg  
EAN of inner box                 4038917084864  
EAN case                         4038917003636

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

