



Product description

Order number	65.724
Product description	Skirted push cap, natural, suitable for tubes Ø 13 mm, 1,000 piece(s)/bag

Product characteristics

Cap	Skirted push cap
Closure type	push cap
Suitable for	tubes Ø 13 mm

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	natural

Purity & certification

Product category	no medical device IVD
Latex-free	yes
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Skirted push cap, natural, suitable for tubes Ø 13 mm

Page 2

Packaging

Minimum order qty.	10000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	10000
Piece(s) / pallet	200000
Depth of case	522 mm
Width of case	322 mm
Height of case	345 mm
Case volume	0.058 cbm
Weight of product	0.001 kg
Weight of case	11.01 kg
EAN of inner box	4038917083720
EAN case	4038917358064

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT