



### **Product description**

Order number	65.790
Product description	Push cap, natural, suitable for tubes Ø 23.5 mm, 1,000 piece(s)/bag

### **Product characteristics**

Cap	standard cap
Closure type	push cap
Suitable for	tubes Ø 23.5 mm

### **Size**

Diameter	22.85 mm
Height of product	22.1 mm

### **Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	natural

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

### **Packaging**

Minimum order qty.                      2000  
Type of smallest subpackaging                      bag  
Piece(s) / inner box                      1000  
Piece (s) / outer case                      2000  
Piece(s) / pallet                          84000  
Depth of case                              388 mm  
Width of case                                378 mm  
Height of case                              250 mm  
Case volume                                0.0367 cbm  
Weight of product                        0.0018 kg  
Weight of case                              4.07 kg  
EAN of inner box                         4038917084130  
EAN case                                    4038917002981

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

