



### **Product description**

Order number 65.793.300  
Product description Push cap, natural, suitable for tubes Ø 16-17 mm, 1,000 piece(s)/bag

### **Product characteristics**

Closure type push cap  
Suitable for tubes Ø 16-17 mm

### **Material & colours**

Product material Low Density Polyethylene (LD-PE)  
Colour of product natural

### **Purity & certification**

Product category no medical device | IVD  
Batched yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	5000
Piece(s) / pallet	200000
Case volume	0.0414 cbm
Weight of product	0.0014 kg
Weight of case	7.5 kg
EAN of inner box	4038917368391
EAN case	4038917368384

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

