



### **Product description**

Order number	65.793.314
Product description	Push cap, yellow, suitable for tubes Ø 16-17 mm, 1,000 piece(s)/bag

### **Product characteristics**

Closure type	push cap
Suitable for	tubes Ø 16-17 mm

### **Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	yellow

### **Purity & certification**

Product category	no medical device   IVD
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	5000
Piece(s) / pallet	200000
Case volume	0.0414 cbm
Weight of product	0.0014 kg
Weight of case	7.58 kg
EAN of inner box	4038917368452
EAN case	4038917368445

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

