



Product description

Order number	65.803.002
Product description	Push cap, green, suitable for tubes Ø 15.7 mm, 1,000 piece(s)/bag

Product characteristics

Cap	standard cap
Closure type	push cap
Suitable for	tubes Ø 15.7 mm

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	green

Purity & certification

Product category	no medical device IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	5000
Piece(s) / pallet	210000
Depth of case	398 mm
Width of case	398 mm
Height of case	250 mm
Case volume	0.0396 cbm
Weight of product	0.0013 kg
Weight of case	7.17 kg
EAN of inner box	4038917121057
EAN case	4038917017626

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

