



Product description

Order number 65.809.305
Product description Push cap, white, suitable for tubes Ø 12 mm, 1,000 piece(s)/bag

Product characteristics

Closure type push cap
Suitable for tubes Ø 12 mm

Material & colours

Product material Low Density Polyethylene (LD-PE)
Colour of product white

Purity & certification

Product category no medical device | IVD
Batched yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	10000
Piece(s) / pallet	400000
Case volume	0.0382 cbm
Weight of product	0.0006 kg
Weight of case	6.82 kg
EAN of inner box	4038917369176
EAN case	4038917369169

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

