



### **Product description**

Order number	55.466
Product description	Tube, working volume: 11.5 ml, (LxØ): 100 x 15.7 mm, material: PS, round base, transparent, push cap, 250 piece(s)/StackPack

### **Product characteristics**

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	11,5 ml
Centrifugation max (RCF)	3500 x g

### **Size**

Diameter	15.7 mm
Length of product	100 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material                      Polystyrene (PS)  
Colour of product                      transparent

### **Purity & certification**

Satisfies the requirement              sealed at 0.6 bar/10 min  
Product category                          no medical device | IVD  
Batched                                      no

### **Packaging**

Minimum order qty.                      1000  
Type of smallest                              StackPack  
subpackaging  
Piece(s) / inner box                      250  
Piece (s) / outer case                      1000  
Piece(s) / pallet                              48000  
Depth of box                                  336 mm  
Width of box                                  102 mm  
Height of box                                  174 mm  
Depth of case                                  388 mm  
Width of case                                  378 mm  
Height of case                                  215 mm  
Case volume                                  0.0315 cbm  
Weight of product                           0.0055 kg  
Weight of case                                6.1 kg  
EAN of inner box                              4038917081740  
EAN case                                        4038917000857

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

