



### **Product description**

Order number	55.473
Product description	Tube, working volume: 8 ml, (LxØ): 84 x 14.5 mm, material: PS, round base, transparent, push cap, 1,000 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	8 ml
Centrifugation max (RCF)	4500 x g

### **Size**

Diameter	14.5 mm
Length of product	84 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

### **Purity & certification**

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	2000
Piece(s) / pallet	60000
Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0.0539 cbm
Weight of product	0.0042 kg
Weight of case	9.1 kg
EAN of inner box	4038917081818
EAN case	4038917000925

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

