



Product description

Order number	55.475.200
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 13 mm, material: PS, round base, transparent, push cap, suitable for Tecan Genesis, FE500, 490 piece(s)/StackPack

Product characteristics

Graduation	no
Closure type	push cap
Suitable for	Tecan Genesis, FE500
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	4000 x g

Size

Diameter	13 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Satisfies the requirement sealed at 0.6 bar/10 min
Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1960
Type of smallest StackPack
subpackaging
Piece(s) / inner box 490
Piece (s) / outer case 1960
Piece(s) / pallet 109760
Depth of box 405 mm
Width of box 77 mm
Height of box 193 mm
Depth of case 428 mm
Width of case 383 mm
Height of case 175 mm
Case volume 0.0287 cbm
Weight of product 0.003 kg
Weight of case 6.83 kg
EAN of inner box 4038917158466
EAN case 4038917158459

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

