



Product description

Order number	55.495
Product description	Tube, working volume: 8.5 ml, (LxØ): 75 x 15.7 mm, material: PS, round base, transparent, push cap, without cap, 1,000 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	8,5 ml
Centrifugation max (RCF)	4000 x g

Size

Diameter	15.7 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Satisfies the requirement sealed at 0.6 bar/10 min
Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 2000
Type of smallest
subpackaging bag
Piece(s) / inner box 1000
Piece (s) / outer case 2000
Piece(s) / pallet 60000
Depth of case 398 mm
Width of case 398 mm
Height of case 340 mm
Case volume 0.0539 cbm
Weight of product 0.0042 kg
Weight of case 9.01 kg
EAN of inner box 4038917082013
EAN case 4038917001120

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

