



### **Product description**

Order number	55.514
Product description	Tube, working volume: 20 ml, (LxØ): 100 x 21.5 mm, material: PP, round base, transparent, push cap, 500 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	20 ml
Centrifugation max (RCF)	4000 x g

### **Size**

Diameter	21.5 mm
Length of product	100 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

### **Purity & certification**

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	24000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0.0657 cbm
Weight of product	0.0067 kg
Weight of case	7.39 kg
EAN of inner box	4038917082099
EAN case	4038917001182

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

