



Product description

Order number	55.516
Product description	Tube, working volume: 8 ml, (LxØ): 100 x 13 mm, material: PP, round base, transparent, push cap, 500 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	8 ml
Centrifugation max (RCF)	3500 x g

Size

Diameter	13 mm
Length of product	100 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Satisfies the requirement sealed at 0.6 bar/10 min
Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 2000
Type of smallest subpackaging bag
Piece(s) / inner box 500
Piece (s) / outer case 2000
Piece(s) / pallet 48000
Depth of case 398 mm
Width of case 398 mm
Height of case 415 mm
Case volume 0.0657 cbm
Weight of product 0.0035 kg
Weight of case 7.79 kg
EAN of inner box 4038917082112
EAN case 4038917001205

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

