



Product description

Order number	55.525.302
Product description	Carrier tube, working volume: 5 ml, (LxØ): 75 x 13 mm, material: PP, round base, brown, push cap, 250 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	5000 x g
Autoclavable (max. °C)	121 °C

Size

Diameter	13 mm
Height of product	75 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product brown

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 250
Type of smallest subpackaging bag
Piece(s) / inner box 250
Piece (s) / outer case 2000
Piece(s) / pallet 60000
Case volume 0.046 cbm
Weight of product 0.0026 kg
Weight of case 5.72 kg
EAN of inner box 4038917367257
EAN case 4038917367240

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

