



Product description

Order number	55.525.505
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 13 mm, material: PP, round base, transparent, push cap, 1,000 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	5000 x g

Size

Diameter	13 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Satisfies the requirement sealed at 0.6 bar/10 min
Product category In vitrodiagnostic device, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 2000
Type of smallest subpackaging bag
Piece(s) / inner box 1000
Piece (s) / outer case 2000
Piece(s) / pallet 84000
Depth of case 398 mm
Width of case 398 mm
Height of case 250 mm
Case volume 0.0396 cbm
Weight of product 0.0026 kg
Weight of case 5.78 kg
EAN of inner box 4038917375955
EAN case 4038917375948

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

