



Product description

Order number	55.526
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PP, round base, transparent, push cap, 500 piece(s)/StackPack

Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	5000 x g

Size

Diameter	12 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Satisfies the requirement sealed at 0.6 bar/10 min
Product category In vitrodiagnostic device, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 4000
Type of smallest StackPack
subpackaging
Piece(s) / inner box 500
Piece (s) / outer case 4000
Piece(s) / pallet 120000
Depth of box 381 mm
Width of box 77 mm
Height of box 176 mm
Depth of case 398 mm
Width of case 358 mm
Height of case 325 mm
Case volume 0.0463 cbm
Weight of product 0.0023 kg
Weight of case 10.17 kg
EAN of inner box 4038917082211
EAN case 4038917001304

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

