



### **Product description**

Order number	55.528
Product description	Tube, working volume: 3 ml, (LxØ): 75 x 10 mm, material: PP, round base, transparent, push cap, 1,000 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	3 ml
Centrifugation max (RCF)	5000 x g

### **Size**

Diameter	10 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent

### **Purity & certification**

Satisfies the requirement sealed at 0.6 bar/10 min  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Batched yes

### **Packaging**

Minimum order qty. 2000  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 1000  
Piece (s) / outer case 2000  
Piece(s) / pallet 96000  
Depth of case 388 mm  
Width of case 378 mm  
Height of case 215 mm  
Case volume 0.0315 cbm  
Weight of product 0.0017 kg  
Weight of case 3.89 kg  
EAN of inner box 4038917082235  
EAN case 4038917001328

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

