



Product description

Order number	55.532.300
Product description	Tube, working volume: 4 ml, (LxØ): 75 x 11.5 mm, material: PP, round base, transparent, push cap, 1,000 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	4 ml
Centrifugation max (RCF)	5000 x g
Autoclavable (max. °C)	121 °C

Size

Diameter	11.5 mm
Height of product	75 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest subpackaging bag
Piece(s) / inner box 1000
Piece (s) / outer case 2000
Piece(s) / pallet 80000
Case volume 0.0382 cbm
Weight of product 0.002 kg
Weight of case 4.77 kg
EAN of inner box 4038917367349
EAN case 4038917367332

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

