



### **Product description**

|                     |                                                                                                                        |
|---------------------|------------------------------------------------------------------------------------------------------------------------|
| Order number        | 55.598.006                                                                                                             |
| Product description | Tube, working volume: 8.5 ml, (LxØ): 75 x 15.7 mm, material: PP, round base, transparent, push cap, 1,000 piece(s)/bag |

### **Product characteristics**

|                          |            |
|--------------------------|------------|
| Graduation               | no         |
| Closure type             | push cap   |
| Base shape               | round base |
| Volume of work           | 8,5 ml     |
| Centrifugation max (RCF) | 4000 x g   |
| Autoclavable (max. °C)   | 121 °C     |

### **Size**

|                   |         |
|-------------------|---------|
| Diameter          | 15.7 mm |
| Height of product | 75 mm   |
| Length of product | 75 mm   |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent

### **Purity & certification**

Product category no medical device | IVD  
Batched yes

### **Packaging**

Minimum order qty. 1000  
Type of smallest subpackaging bag  
Piece(s) / inner box 1000  
Piece (s) / outer case 2000  
Piece(s) / pallet 56000  
Case volume 0.0649 cbm  
Weight of product 0.0024 kg  
Weight of case 5.9 kg  
EAN of inner box 4038917363136  
EAN case 4038917363129

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

