



Product description

Order number 58.485
Product description Tube, working volume: 7 ml, (LxØ): 50 x 16 mm, material: PS, flat base, transparent, push cap, 1,000 piece(s)/bag

Product characteristics

Graduation no
Closure type push cap
Base shape flat base
Volume of work 7 ml
Centrifugation max (RCF) 2500 x g

Size

Diameter 16 mm
Length of product 50 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	2000
Piece(s) / pallet	72000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0028 kg
Weight of case	6.34 kg
EAN of inner box	4038917081924
EAN case	4038917001038

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

