



Product description

| | |
|---------------------|--|
| Order number | 58.618 |
| Product description | Tube, working volume: 12 ml, (LxØ): 95 x 16.5 mm, material: PS, skirted conical base, transparent, push cap, without cap, 500 piece(s)/bag |

Product characteristics

| | |
|----------------|----------------------|
| Graduation | no |
| Closure type | push cap |
| Base shape | skirted conical base |
| Volume of work | 12 ml |

Size

| | |
|-------------------|---------|
| Diameter | 16.5 mm |
| Length of product | 95 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 1000000
Type of smallest
subpackaging bag
Piece(s) / inner box 500
Piece (s) / outer case 1000
Piece(s) / pallet 36000
Case volume 0.0492 cbm
Weight of product 0.0005 kg
EAN of inner box 4038917082938
EAN case 4038917001922

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

