



### **Product description**

Order number 58.631  
Product description Tube, working volume: 23 ml, (LxØ): 75 x 23.5 mm, material: PP, flat base, transparent, push cap, 500 piece(s)/bag

### **Product characteristics**

Graduation no  
Closure type push cap  
Base shape flat base  
Volume of work 23 ml  
Centrifugation max (RCF) 3500 x g

### **Size**

Diameter 23.5 mm  
Length of product 75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

### **Purity & certification**

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	24000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0.0657 cbm
Weight of product	0.0055 kg
Weight of case	6.18 kg
EAN of inner box	4038917083027
EAN case	4038917002004

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

