



## Product Specification

**Screw cap tube, 5 ml, (LxØ): 75 x 13 mm, round base, PP, without cap, 1,000 piece(s)/bag**

**Page 1**



### Product description

|                     |  |
|---------------------|--|
| Order number        | 60.504.010   |
| Product description | Screw cap tube, working volume: 5 ml, (LxØ): 75 x 13 mm, round base, transparent, material: PP, without cap, 1,000 piece(s)/bag, 2,000 piece(s)/carton |

### Product characteristics

|                          |            |
|--------------------------|------------|
| Graduation               | no         |
| Closure type             | screw cap  |
| Base shape               | round base |
| Volume of work           | 5 ml       |
| Centrifugation max (RCF) | 5000 x g   |

### Size

|                   |       |
|-------------------|-------|
| Diameter          | 13 mm |
| Length of product | 75 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Screw cap tube, 5 ml, (LxØ): 75 x 13 mm, round base, PP, without cap, 1,000 piece(s)/bag**

**Page 2**

### Material & colours

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | transparent        |
| Tube material     | Polypropylene (PP) |

### Purity & certification

|                           |                               |
|---------------------------|-------------------------------|
| Satisfies the requirement | IATA                          |
| Product category          | In vitrodiagnostic device, CE |
| CE certified              | CE                            |
| Batched                   | yes                           |

### Packaging

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 2000          |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 1000          |
| Piece (s) / outer case        | 2000          |
| Piece(s) / pallet             | 72000         |
| Depth of case                 | 398 mm        |
| Width of case                 | 398 mm        |
| Height of case                | 265 mm        |
| Case volume                   | 0.042 cbm     |
| Weight of product             | 0.0024 kg     |
| Weight of case                | 5.43 kg       |
| EAN of inner box              | 4038917309059 |
| EAN case                      | 4038917531764 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**