



**Product Specification**  
**Screw cap tube, 5 ml, (LxØ): 75 x 13 mm, round base, PP, without cap, 1,000 piece(s)/bag**

**Page 1**



**Product description**

Order number	60.504.021
Product description	Screw cap tube, working volume: 5 ml, (LxØ): 75 x 13 mm, round base, transparent, material: PP, without cap, 1,000 piece(s)/bag, 2,000 piece(s)/carton

**Product characteristics**

Graduation	no
Closure type	screw cap
Base shape	round base
Volume of work	5 ml

**Size**

Diameter	13 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap tube, 5 ml, (LxØ): 75 x 13 mm, round base, PP, without cap, 1,000 piece(s)/bag**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Tube material	Polypropylene (PP)

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	2000
Piece(s) / pallet	72000
Case volume	0.0492 cbm
Weight of product	0.0024 kg
Weight of case	5.74 kg
EAN of inner box	4038917348201
EAN case	4038917348195

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

