



### **Product description**

Order number	60.540.500
Product description	Screw cap tube, working volume: 13 ml, (LxØ): 101 x 16.5 mm, material: PP, round base, transparent, screw cap, natural, cap assembled, sterile, 500 piece(s)/bag

### **Product characteristics**

Graduation	no
Cap	standard screw cap
Closure type	screw cap
Base shape	round base
Volume of work	13 ml
Centrifugation max (RCF)	4000 x g

### **Size**

Diameter	16.5 mm
Length of product	101 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural

### **Purity & certification**

Satisfies the requirement	IATA
Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### **Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	500
Piece(s) / pallet	28000
Depth of case	598 mm
Width of case	398 mm
Height of case	125 mm
Case volume	0.0298 cbm
Weight of product	0.0059 kg
Weight of case	3.52 kg
EAN of inner box	4038917110495
EAN case	4038917012881

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

