



### Product description

|                     |  |
|---------------------|--|
| Order number        | 60.541.009   |
| Product description | Screw cap tube, working volume: 13 ml, (LxØ): 101 x 16.5 mm, material: PP, round base with skirted base, transparent, screw cap, without cap, 500 piece(s)/bag |

### Product characteristics

|                          |                              |
|--------------------------|------------------------------|
| Graduation               | no                           |
| Closure type             | screw cap                    |
| Base shape               | round base with skirted base |
| Volume of work           | 13 ml                        |
| Centrifugation max (RCF) | 4000 x g                     |
| Autoclavable (max. °C)   | 121 °C                       |

### Size

|                   |         |
|-------------------|---------|
| Diameter          | 16.5 mm |
| Length of product | 101 mm  |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | transparent        |
| Tube material     | Polypropylene (PP) |

### **Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Batched          | yes                     |

### **Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1000          |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 500           |
| Piece (s) / outer case        | 1000          |
| Piece(s) / pallet             | 32000         |
| Case volume                   | 0.0558 cbm    |
| Weight of product             | 0.0067 kg     |
| Weight of case                | 7.2 kg        |
| EAN of inner box              | 4038917110563 |
| EAN case                      | 4038917032438 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

