



### **Product description**

Order number	60.542
Product description	Screw cap tube, working volume: 8 ml, (LxØ): 57 x 16.5 mm, material: PP, flat base, transparent, screw cap, natural, cap enclosed, 1,000 piece(s)/bag

### **Product characteristics**

Graduation	no
Cap	standard screw cap
Closure type	screw cap
Base shape	flat base
Volume of work	8 ml
Centrifugation max (RCF)	4000 x g

### **Size**

Diameter	16.5 mm
Length of product	57 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural

### **Purity & certification**

Satisfies the requirement	IATA
Product category	no medical device   IVD
Batched	yes

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	1000
Piece(s) / pallet	42000
Depth of case	388 mm
Width of case	378 mm
Height of case	250 mm
Case volume	0.0367 cbm
Weight of product	0.0042 kg
Weight of case	4.64 kg
EAN of inner box	4038917082358
EAN case	4038917001458

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

