



### **Product description**

Order number	60.542.324
Product description	Screw cap tube, working volume: 8 ml, (LxØ): 57 x 16.5 mm, material: PP, flat base, transparent, screw cap, natural, cap assembled, sterile, 100 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	screw cap
Base shape	flat base
Volume of work	8 ml
Centrifugation max (RCF)	4000 x g

### **Size**

Diameter	16.5 mm
Length of product	57 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural

### **Purity & certification**

Satisfies the requirement	IATA
Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	36000
Case volume	0.0436 cbm
Weight of product	0.0044 kg
Weight of case	4.79 kg
EAN of inner box	4038917367851
EAN case	4038917367844

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

