



### **Product description**

Order number	60.551.023
Product description	Screw cap tube, working volume: 10 ml, (LxØ): 79 x 16 mm, material: PP, round base with skirted base, transparent, screw cap, without cap, 500 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	screw cap
Base shape	round base with skirted base
Volume of work	10 ml
Centrifugation max (RCF)	4000 x g
Autoclavable (max. °C)	121 °C

### **Size**

Diameter	16 mm
Height of product	79 mm
Length of product	79 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent

### **Purity & certification**

Product category no medical device | IVD  
Batched yes

### **Packaging**

Minimum order qty. 2000  
Type of smallest subpackaging bag  
Piece(s) / inner box 500  
Piece (s) / outer case 2000  
Piece(s) / pallet 32000  
Case volume 0.0813 cbm  
Weight of product 0.0024 kg  
Weight of case 5.9 kg  
EAN of inner box 4038917317801  
EAN case 4038917317795

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

