



### **Product description**

Order number	60.732
Product description	Screw cap tube, working volume: 15 ml, (LxØ): 76 x 20 mm, material: PP, skirted conical base, transparent, screw cap, natural, cap enclosed, 100 piece(s)/bag

### **Product characteristics**

Graduation	no
Cap	standard screw cap
Closure type	screw cap
Base shape	skirted conical base
Volume of work	15 ml
Centrifugation max (RCF)	4500 x g

### **Size**

Diameter	20 mm
Length of product	76 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural

### **Purity & certification**

Satisfies the requirement	ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

### **Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	500
Piece(s) / pallet	18000
Depth of case	398 mm
Width of case	398 mm
Height of case	265 mm
Case volume	0.042 cbm
Weight of product	0.0064 kg
Weight of case	3.76 kg
EAN of inner box	4038917083805
EAN case	4038917002691

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

