



### **Product description**

Order number	60.9921.525
Product description	Screw cap tube, working volume: 5 ml, (LxØ): 50 x 16 mm, material: PS, flat base, transparent, 250 piece(s)/bag

### **Product characteristics**

Graduation	no
Base shape	flat base
Volume of work	5 ml
Centrifugation max (RCF)	5000 x g
Centrifugation temperature	21
Centrifugation time	15 Min

### **Size**

Diameter	16 mm
Length of product	50 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent
Tube material	Polystyrene (PS)

### **Purity & certification**

Satisfies the requirement	sealed at 18"hg/10min
Latex-free	yes
Pyrogen-free	no
Safety data sheet	no

### **Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	2000
Piece(s) / pallet	96000
Depth of case	608 mm
Width of case	264 mm
Height of case	324 mm
Case volume	0.0519 cbm
Weight of case	7 kg
EAN of inner box	4038917150811
EAN case	4038917037167

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

