



Product description

Order number	62.467.012
Product description	Tube, working volume: 8 ml, (LxØ): 100 x 13 mm, material: PS, round base, transparent, push cap, with print, label/print: white, with graduations, 500 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	push cap
Base shape	round base
Volume of work	8 ml
Centrifugation max (RCF)	3500 x g

Size

Diameter	13 mm
Height of product	100 mm
Length of product	100 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 2000
Type of smallest
subpackaging bag
Piece(s) / inner box 500
Piece (s) / outer case 2000
Piece(s) / pallet 56000
Case volume 0.0649 cbm
Weight of product 0.0041 kg
Weight of case 9.2 kg
EAN of inner box 4038917106856
EAN case 4038917032599

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

