




Product Specification
Screw cap tube, 5 ml, (LxØ): 75 x 13 mm, round base, PP, cap enclosed, 1,000 piece(s)/bag



Product description

| | |
|---------------------|---|
| Order number | 62.504.016 |
| Product description | Screw cap tube, working volume: 5 ml, (LxØ): 75 x 13 mm, round base, transparent, material: PP, with print, label/print: white, with graduations, cap enclosed, 1,000 piece(s)/bag, 1,000 piece(s)/carton |

Product characteristics

| | |
|-----------------------|------------|
| Label/ Print | with print |
| Colour of print/label | white |
| Graduation | yes |
| Closure type | screw cap |
| Base shape | round base |
| Volume of work | 5 ml |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

Screw cap tube, 5 ml, (LxØ): 75 x 13 mm, round base, PP, cap enclosed, 1,000 piece(s)/bag

Page 2

Size

| | |
|-------------------|-------|
| Diameter | 13 mm |
| Length of product | 75 mm |

Material & colours

| | |
|-------------------|-----------------------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Tube material | Polypropylene (PP) |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 1000 |
| Piece (s) / outer case | 1000 |
| Piece(s) / pallet | 56000 |
| Case volume | 0.029 cbm |
| Weight of product | 0.0039 kg |
| Weight of case | 4.38 kg |
| EAN of inner box | 4038917348058 |
| EAN case | 4038917348041 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT