



**Product Specification**  
**Screw cap tube, 5 ml, (LxØ): 75 x 13 mm, round base, PP, cap enclosed, 100 piece(s)/bag**

**Page 1**



**Product description**

Order number	62.504.039
Product description	Screw cap tube, working volume: 5 ml, (LxØ): 75 x 13 mm, round base, brown, material: PP, with print, label/print: black, with graduations, cap enclosed, brown, 100 piece(s)/bag, 1,000 piece(s)/carton

**Product characteristics**

Label/ Print	with print
Colour of print/label	black
Graduation	yes
Closure type	screw cap
Light protection	yes
Base shape	round base
Volume of work	5 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

Screw cap tube, 5 ml, (LxØ): 75 x 13 mm, round base, PP, cap enclosed, 100 piece(s)/bag

Page 2

### Size

Diameter	13 mm
Length of product	75 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	brown
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	brown
Tube material	Polypropylene (PP)

### Purity & certification

Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

### Packaging

Minimum order qty.	50000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	40000
Depth of case	458 mm
Width of case	408 mm
Height of case	205 mm
Case volume	0.0383 cbm
Weight of product	0.0039 kg
Weight of case	4.39 kg
EAN of inner box	4038917413008
EAN case	4038917412995

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**