



Product description

Order number	62.536.009
Product description	Tube, working volume: 7 ml, (LxØ): 50 x 16 mm, material: PP, flat base, transparent, push cap, with print, label/print: white, with graduations, 1,000 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	push cap
Base shape	flat base
Volume of work	7 ml
Autoclavable (max. °C)	121 °C

Size

Diameter	16 mm
Height of product	50 mm
Length of product	50 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Satisfies the requirement IATA
Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest subpackaging bag
Piece(s) / inner box 1000
Piece (s) / outer case 2000
Piece(s) / pallet 72000
Case volume 0.0492 cbm
Weight of product 0.002 kg
Weight of case 5.1 kg
EAN of inner box 4038917110037
EAN case 4038917033107

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

