



**Product Specification**  
**Screw cap tube, 30 ml, (LxØ): 107 x 25 mm, PP,  
with print**



**Product description**

Order number	62.543.001
Product description	Screw cap tube, working volume: 30 ml, (LxØ): 107 x 25 mm, material: PP, skirted conical base, transparent, screw cap, natural, cap assembled, with print, label/print: black, with graduations, sterile, 50 piece(s)/bag

**Product characteristics**

Label/ Print	with print
Colour of print/label	black
Graduation	yes
Cap	standard screw cap
Closure type	screw cap
Base shape	skirted conical base
Volume of work	30 ml
Centrifugation max (RCF)	8800 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap tube, 30 ml, (LxØ): 107 x 25 mm, PP,  
with print**

**Size**

Diameter	25 mm
Length of product	107 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural

**Purity & certification**

Satisfies the requirement	ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Minimum order qty.	400
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	400
Piece(s) / pallet	11200
Depth of case	588 mm
Width of case	398 mm
Height of case	255 mm
Case volume	0.0597 cbm
Weight of product	0.0122 kg
Weight of case	5.59 kg
EAN of inner box	4038917111638
EAN case	4038917013284

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

