



**Product Specification**  
**Screw cap tube, 30 ml, (LxØ): 107 x 25 mm, PP,  
with print**



**Product description**

Order number	62.543.030
Product description	Screw cap tube, working volume: 30 ml, (LxØ): 107 x 25 mm, material: PP, skirted conical base, white, screw cap, natural, cap assembled, with print, label/print: black, with graduations, 125 piece(s)/bag

**Product characteristics**

Label/ Print	with print
Colour of print/label	black
Graduation	yes
Closure type	screw cap
Light protection	yes
Base shape	skirted conical base
Volume of work	30 ml
Centrifugation max (RCF)	8800 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap tube, 30 ml, (LxØ): 107 x 25 mm, PP,  
with print**

**Page 2**

**Size**

Diameter	25 mm
Length of product	107 mm
Length including cap	110 mm
Length excluding cap	107 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	white
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural

**Purity & certification**

Satisfies the requirement	ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	125
Piece (s) / outer case	500
Piece(s) / pallet	12000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0.0657 cbm
Weight of product	0.0122 kg
Weight of case	6.84 kg
EAN of inner box	4038917111713
EAN case	4038917076203

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**