



Product Specification
**Screw cap tube, 50 ml, (LxØ): 114 x 28 mm, PP,
with print**

Page 1



Product description

Order number	62.547.275
Product description	Screw cap tube, working volume: 50 ml, (LxØ): 114 x 28 mm, material: PP, conical base, transparent, HD-PE, cap: natural, cap assembled, with print, label/print: white/blue, with graduations, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 25 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white/blue
Graduation	yes
Closure type	screw cap
Base shape	conical base
Volume of work	50 ml
Centrifugation max (RCF)	20000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
**Screw cap tube, 50 ml, (LxØ): 114 x 28 mm, PP,
with print**

Size

Diameter	28 mm
Length of product	114 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural

Purity & certification

Satisfies the requirement	IATA
Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	250
Type of smallest subpackaging	bag
Piece(s) / inner box	25
Piece (s) / outer case	250
Piece(s) / pallet	9000
Case volume	0.0492 cbm
Weight of product	0.0161 kg
Weight of case	4.15 kg
EAN of inner box	4038917387330
EAN case	4038917387323

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

