



Product Specification
**Screw cap tube, 10 ml, (LxØ): 79 x 16 mm, PP,
with print**



Product description

Order number	62.551.201
Product description	Screw cap tube, working volume: 10 ml, (LxØ): 79 x 16 mm, material: PP, round base with skirted base, transparent, screw cap, natural, cap assembled, with print, label/print: white, with graduations, sterile, 100 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	screw cap
Base shape	round base with skirted base
Volume of work	10 ml
Centrifugation max (RCF)	4000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Screw cap tube, 10 ml, (LxØ): 79 x 16 mm, PP,
with print**

Page 2

Size

Diameter 16 mm
Length of product 79 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap natural

Purity & certification

Satisfies the requirement IATA
Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest subpackaging bag
Piece(s) / inner box 100
Piece (s) / outer case 1000
Piece(s) / pallet 36000
Depth of case 413 mm
Width of case 257 mm
Height of case 464 mm
Case volume 0.0492 cbm
Weight of product 0.0044 kg
Weight of case 5.04 kg
EAN of inner box 4038917112901
EAN case 4038917072113

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT