




Product Specification
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PS,
with print**



Product description

Order number	62.553.001
Product description	Screw cap tube, working volume: 15 ml, (LxØ): 120 x 17 mm, material: PS, conical base, transparent, screw cap, red, cap enclosed, with print, label/print: white/blue, with graduations, 50 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white/blue
Graduation	yes
Closure type	screw cap
Base shape	conical base
Volume of work	15 ml
Centrifugation max (RCF)	4100 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PS,
with print**

Size

Diameter 17 mm
Length of product 120 mm

Material & colours

Product material Polystyrene (PS)
Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap red

Purity & certification

Satisfies the requirement IATA
Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 500
Type of smallest subpackaging bag
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 20000
Case volume 0.0382 cbm
Weight of product 0.0086 kg
Weight of case 4.54 kg
EAN of inner box 4038917112932
EAN case 4038917033428

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

