



Product Specification
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PS,
with print**



Product description

Order number	62.553.042
Product description	Screw cap tube, working volume: 15 ml, (LxØ): 120 x 17 mm, material: PS, conical base, transparent, screw cap, yellow, cap assembled, with print, label/print: white/blue, with graduations, sterile, 50 piece(s)/polystyrene container

Product characteristics

Label/ Print	with print
Colour of print/label	white/blue
Graduation	yes
Closure type	screw cap
Base shape	conical base
Volume of work	15 ml
Centrifugation max (RCF)	4100 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PS,
with print**

Size

Diameter	17 mm
Length of product	120 mm

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	polystyrene container
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	16000
Depth of case	648 mm
Width of case	420 mm
Height of case	145 mm
Case volume	0.0395 cbm
Weight of product	0.0074 kg
Weight of case	4.64 kg
EAN of inner box	4038917112970
EAN case	4038917013864

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

