



**Product Specification**  
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PP,  
with print**



**Product description**

Order number	62.554.002
Product description	Screw cap tube, working volume: 15 ml, (LxØ): 120 x 17 mm, material: PP, conical base, transparent, screw cap, red, cap assembled, with print, label/print: white/blue, with graduations, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 50 piece(s)/polystyrene container

**Product characteristics**

Label/ Print	with print
Colour of print/label	white/blue
Graduation	yes
Cap	standard screw cap
Closure type	screw cap
Base shape	conical base
Volume of work	15 ml
Centrifugation max (RCF)	20000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PP,  
with print**

**Size**

Diameter 17 mm  
Length of product 120 mm

**Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap red

**Purity & certification**

Satisfies the requirement IATA  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Purity standard free from DNA/DNase/RNase, free from  
pyrogens/endotoxins, non-cytotoxic, sterile  
Sterilisation Electron irradiation  
Pyrogen-free yes  
Batched yes

**Packaging**

Minimum order qty. 500  
Type of smallest polystyrene container  
subpackaging  
Piece(s) / inner box 50  
Piece (s) / outer case 500  
Piece(s) / pallet 16000  
Depth of case 435 mm  
Width of case 143 mm  
Height of case 702 mm  
Case volume 0.0436 cbm  
Weight of product 0.0057 kg  
Weight of case 4.13 kg  
EAN of inner box 4038917565103  
EAN case 4038917565110

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

