




Product Specification
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PP,
with print**

Page 1



Product description

Order number	62.554.016
Product description	Screw cap tube, working volume: 15 ml, (LxØ): 120 x 17 mm, material: PP, conical base, transparent, screw cap, natural, cap assembled, with print, label/print: white/blue, with graduations, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 50 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white/blue
Graduation	yes
Closure type	screw cap
Base shape	conical base
Volume of work	15 ml
Centrifugation max (RCF)	20000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PP,
with print**

Size

Diameter 17 mm
Length of product 120 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap natural

Purity & certification

Satisfies the requirement IATA
Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard free from DNA/DNase/RNase, free from
pyrogens/endotoxins, non-cytotoxic, sterile
Sterilisation Electron irradiation
Pyrogen-free yes
Batched yes

Packaging

Minimum order qty. 500
Type of smallest bag
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 20000
Depth of case 598 mm
Width of case 398 mm
Height of case 165 mm
Case volume 0.0393 cbm
Weight of product 0.0063 kg
Weight of case 3.74 kg
EAN of inner box 4038917233538
EAN case 4038917233521

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

