



Product Specification
**Screw cap tube, 30 ml, (LxØ): 84 x 30 mm, PP,
with print**

Page 1



Product description

Order number	62.555.001
Product description	Screw cap tube, working volume: 30 ml, (LxØ): 84 x 30 mm, material: PP, flat base, transparent, screw cap, natural, cap assembled, with print, label/print: white, with graduations, sterile, 50 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	screw cap
Base shape	flat base
Volume of work	30 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Screw cap tube, 30 ml, (LxØ): 84 x 30 mm, PP,
with print**

Page 2

Size

Diameter	30 mm
Length of product	84 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	natural

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	450
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	450
Piece(s) / pallet	12600
Depth of case	588 mm
Width of case	398 mm
Height of case	255 mm
Case volume	0.0597 cbm
Weight of product	0.011 kg
Weight of case	5.71 kg
EAN of inner box	4038917113205
EAN case	4038917013949

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT